

Intelligent Design

The core product of the RescueNet ePCR Suite is TabletPCR. Its user interface brings the latest ideas and technology in usability to the electronic field data experience.

- Optimized for touch-screen tablets
- Easily tab from field-to-field when using keyboard-based input
- QuickLog allows fast time-stamping of common procedures
- Complete Call Rules with hyperlinks to required fields
- Spell-check with a medical dictionary to help ensure professional reports.


RescueNet TabletPCR is the most effective field data solution for record accuracy, data integration and advanced analysis. TabletPCR is widely used in the EMS industry and has consistently improved from years of feedback from the field.


- Prevent incomplete patient care reports
- Improve medic efficacy
- Reduce risk in the event of litigation or audits
- Eliminate inefficient paper processes
- Ensure protocol is followed consistently
- Utilize integrated field data for syndromic surveillance.

Configurable to Meet Your Needs

With configurable screens for TabletPCR and WebPCR, you can quickly customize the system to meet the needs of your organization. Administrators can modify screen layout, change labels and create their own pick lists – all without costly custom programming. Complete Call Rules ensure that calls are both accurate and complete.

In simple terms, the software adapts to your business – your business doesn't have to adapt to the software. Designed for easy data entry, training time to transition from paper systems is minimal. The system administrator can customize the data entry screens to add fields, change label names and rearrange fields on the screen (see image).





ePCR Suite Modules and Features

The RescueNet ePCR Suite has a number of features and components that give you control over the documentation process in your organization:

RescueNet WebPCR

- Web-based PCR editor provides PCR editing and review from any workstation with a compatible web browser
- Built with the same screen layout as TabletPCR so that everything is in the same location in the application as it is with TabletPCR
- Integral component of the PCR Workflow module, which allows the automatic routing of PCRs for review and the manual routing of PCRs for corrections

- Allows routing of calls to users on the system for review or correction
- Supports "route and return" to ensure that after the PCR is reviewed, it is returned to sender.

RescueNet eDistribution

- Saves money by allowing the automatic distribution of files and updates to mobile units without having to physically "touch" the mobile units.

RescueNet ePCR Workflow Editor

- Advanced, graphical editor to design the flow of PCRs through the review process
- Designed to allow quick layout of conceptual flow to simplify the design process.

Reporting on Results

RescueNet Reporting provides administrative personnel with out-of-the box reports on many aspects of the operation. You can use these reports to collect and print ePCR information. In addition to the standard reports provided with RescueNet TabletPCR, you can also create custom reports using Crystal Reports' award-winning report designer. Over fifty reports are included in the following categories:

- Administration reports
- Treatments and impressions reports
- Medic skills reports
- Trip lists
- Performance > On Scene reports
- General > HIPAA reports
- Statistics reports
- RescueNet Dispatch and Billing reports



ISV/Software Solutions

www.zolldata.com

US Sales 800 474 4489 | US Fax 303 801 0001 | US Direct Dial +1 303 801 0000
Asia/Pacific +61 2 9420 8733 | Europe +44 (0) 1925 846400